

# Foot Buckets and Calf/Legrests

## 1 Calf/Legrest Pad

The calf/legrest pads can be attached to most wheelchairs by using a 434H1 S.T.S.

1

## 2 Foot Bucket Cushions

Foot bucket cushions provide comfortable, safe support and are compatible with most wheelchairs on the market.

2






### Calf / Legrest Pad

The calf/legrest pads provide an ideal anatomic fit and are cushioned for maximum comfort. The pads attach to most wheelchairs by using a 434H1 S.T.S. clamp and 434F3 rod assembly.

Article no.	Size
440P2=1-7	Width 6 11/16" (17cm) Height 9 1/16" (23cm)



### Dropped S.T.S. Mounting Kit

Article no.	Color	Side
434F3=L	Chrome	Left
434F3=R	Chrome	Right
434F3=L-7	Black	Left
434F3=R-7	Black	Right

Size A	Size B	Size C
4 5/16" (11cm)	3 3/4" (9.5cm)	1 3/4" (4.5cm)



### S.T.S. Clamps

These clamps can be fastened to the wheelchair without alterations. Please state tubing size when ordering. For accuracy, please contact wheelchair manufacturer for tubing size. For Ottobock wheelchairs, you need S.T.S. clamps for tube diameter 1" (2.5cm).

Article no.	Color	Tube Ø
434H1=3/4	Chrome	3/4" (1.9cm)
434H1=7/8	Chrome	7/8" (2.2cm)
434H1=1	Chrome	1" (2.5cm)
434H1=3/4-7	Black	3/4" (1.9cm)
434H1=7/8-7	Black	7/8" (2.2cm)
434H1=1-7	Black	1" (2.5cm)



### Foot Bucket Shells

Consisting of an ABS bucket and a soft polyurethane foam padding, foot buckets offer comfortable and secure support. They are mounted directly onto existing footplates of standard wheelchairs and can be used for either seating systems or wheelchairs. Padding and buckets must be ordered separately. Foot buckets come in two sizes (large and medium – both with a 90 degree angle).

Article no.	Width	Height	Depth
440R9=2-7 (medium)	8 11/16" (22cm)	10" (25.5cm)	5 7/8" (15cm)
440R9=1-7 (large)	12 4/5" (32.5cm)	10" (25.5cm)	7 1/8" (18.5cm)



### Foot Bucket Cushions

Foot bucket cushions are used with Ottobock seating systems or standard wheelchairs to provide comfortable and safe support. The cushion is made from a polyurethane foam that is contoured to comfortably position the foot.

Article no.	Width	Height	Depth
440C2=2-7 (medium)	8 11/16" (22cm)	10" (25.5cm)	5 7/8" (15cm)
440C2=1-7 (large)	12 4/5" (32.5cm)	10" (25.5cm)	7 1/8" (18.5cm)



### Abductor

In combination with the mounting kit, the abductor can be flipped down by pressing the push button. It can only be attached to firm seat frames or seat plates.

Article no.	Width	Height	Depth
446C1=1-7	3" (13.5cm)	3 15/16" (10cm)	3 15/16" (10cm)



### Mounting Kit for Abductor

Article no.	Description
438H2	Mounting Kit for Abductor